
## **Terraillon**.

## Webcoach Prime Android

- Download the Terraillon Wellness Coach app from Google Play Store (Android). Create your account (if you don't already have one). Ensure your Bluetooth connection is on.
- **2.** Hold the UNIT button down at the back of the WEB COACH PRIME scale until the display shows a rotating symbol to start pair up.
- **3.** Access "SYNCHRONISE MY DEVICES" (fig.1) from the menu on the right hand corner. Tap "START SYNCRONISATION" (fig. 2) from your phone (this may take a few seconds). Once the operation succeeds and the Web Coach symbol appears, tap "SYNCHRONISE ALL DEVICES" (fig.3).



**4.** The 2 devices should now be synchronised. (fig.4). It is important to bind your account with the WEB COACH PRIME scale for automatic recognition of user during weighing. Tap "NEXT". Choose a profile and edit the profile with your name (fig.5). Then press "CONFIRM".

米 👀 宮 📶 44	4% 🖬 13:47	**	¥ 🖸 😤 📶 44%
$\mathbf{k}$ Terraillon. $\equiv$		Bind a user account	
	s	Please select a profile to synchr	onize your scale
		P1 : Profile nº1	
		P2 : Profile nº2	
	Ð	P3: Profile nº3	
		P4 : Profile nº4	
SYNCHR, SE ALL DEVICES		P5 : Profile n°5	
Synchronising devices Device 1/1: Web Coach		P6 : Profile n°6	
		P7 : Profile n°7	
		P8 : Profile nº8	
		Cancel	



## **Terraillon**®

- **5.** Once you have bind a user account, your WEB COACH PRIME is set and ready to be used (on left corner you see PX, with X being your profile number). If you are not doing your first scale right after the synchronisation, the display will stay off. Select your profile name again and press "SET" (fig.5).
- **6.** Use the WEB COACH PRIME scale on a hard, flat surface. Step onto the scale and you will see your weight displayed on the screen. Once your weight is stable, the display blinks.
  - If you step out right away, only your weight will be measured.
  - If you stay longer on the scale, the "BIA" symbol appears on the left corner of the display. It means it starts analysing your body composition.

In both cases, a Bluetooth icon appears when weighting is complete.

7. A log book icon will flash on the screen confirming that the connection with the app is complete and that the data has been transmitted to the Wellness Coach app. From the Menu on the right hand corner, select 'ACCESS MY PROGRESS CHART'. Your dashboard will appear. Scroll down until you see the 'MY WEIGHT' and you will see notifications per weigh in session. Review your results and either save or delete connections (fig.6).



Please note that the Wellness Coach APP is supported with Android 4.3 and above.