# *BARBER***BOSS**

## **QR-3641**



Importer: Elite Capital Limited

Address: Studio 4 Clarks Courtyard, 145 Granville Street,  
Birmingham, B1 1SB, United Kingdom.

Website: [barberboss.co.uk](http://barberboss.co.uk)

Email: [care@barberboss.co.uk](mailto:care@barberboss.co.uk)

QR-3641

8

9

10

11


12


13

14

15

Welcome to BarberBoss and congratulations on your new purchase! To ensure your product is fully protected, please activate your 12-month warranty within 90 days by visiting: [www.barberboss.co.uk/barberboss-product-warranty](http://www.barberboss.co.uk/barberboss-product-warranty). Beyond our warranty, you are also covered by your country's legal consumer rights, ensuring complete protection for your purchase.

1. Massage Head
2. LED Display
3. Speed Button
4. Power Button
5. Charging Interface
6. Silicone Handle
7. Charging Cable
8. Flat Head (Suitable for large flat muscle groups such as the back, waist, and arms)
9. Sphere Head (Ideal for large muscle groups across the body)
10. Cylinder Head (For joints and targeted deep-tissue areas)
11. U-Head (For the spine and area behind the ankle)
12. Finger Head (For deep-tissue area such as soles, joints, palms)
13. Silicone Swamm Head (For the whole body, especially medium to large muscle areas, due to its soft and gentle surface)
14. Round Scoop Head (For small muscle areas)
15. Camber Head (Ideal for the legs, arms, and neck)

#### Product Details

|                |                           |
|----------------|---------------------------|
| Product Name   | QR-3641 Small Massage Gun |
| Model Number   | QR-3641                   |
| Charging Input | DC 5V 1A                  |
| Charging Time  | ≤4 Hours                  |
| Usage Time     | 80 mins (Max level)       |

#### Charging Instructions:

- Make sure the device is switched OFF before charging.
- Charge the product fully before you use it for the first time and after a long period of disuse. For optimal performance, it is best to use the device when fully charged.
- When the device power is low, the LED Display will flash to remind you to charge the device.
- When the device is charging, the LED display lights will blink in a form of circle. Once the battery is fully charged, all lights in the circle will stay on continuously.

#### How to use:

- Long press the power button to Turn ON or Turn OFF the device.
- Press '+' button to increase massage force (from 1- 9). The LED Display will correspond to the selected level of force.

- When the massage gun is pressed too hard against the body, the force indicator light will turn red (it will first turn yellow as a warning).

#### **Precautions:**

- Charge only with DC 5V 1A certified adapter.
- Do not overcharge or charge near water. Do not expose to water or moisture. This device is not waterproof.
- Do not use if damaged or malfunctioning.
- The product should be used by adults only and keep out of reach of children.
- For your safety, do not use the massage gun while it is charging.
- Use the device only on soft tissues of the body without causing pain or discomfort. Do not use it on head, bones of sensitive parts of the body.
- Stop using the device once there is any sign of pain or discomfort.
- Only the attachments on the recommended muscle groups.
- To protect the bones and muscles, avoid long-term fixed-point massage.
- Do not immerse the device in water or allow water to enter the device.
- Do not use on injured or swollen areas, or if you have open wounds, bruises, or skin conditions.
- Consult your doctor before use if you have any medical conditions, such as heart problems, diabetes, or if you are pregnant.
- Do not use for more than 60 seconds on a single area. Move the device slowly across muscles.
- Do not operate while sleeping, drowsy, or under the influence of drugs or alcohol.
- Use only the provided charger and accessories.
- Stop use immediately if you feel pain, discomfort, dizziness, or notice any abnormal reaction.
- Do not disassemble or attempt to repair the device yourself. Contact customer service for assistance.

#### **Cleaning:**

- Power off and disconnect before cleaning.
- To clean, wipe the surface of the device with a damp towel and dry it with a soft towel.

#### **Customer Service**

- Detailed Instructions, Cleaning, Troubleshooting, Environmental Protection and Recycling: <https://www.barberboss.co.uk/a/FAQ>
- Replacement parts & accessories: [www.barberboss.co.uk/shop](http://www.barberboss.co.uk/shop)
- Email: [care@barberboss.co.uk](mailto:care@barberboss.co.uk)
- Phone/WhatsApp: 0800 612 6688 (UK), +44800 612 6688 (International)
- Tutorials: [youtube.com/@BarberBossUK](https://youtube.com/@BarberBossUK)

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