*BARBER***BOSS**

QR-7116


Importer: Elite Capital Limited

Address: Studio 4 Clarks Courtyard, 145 Granville Street,
Birmingham, B1 1SB, United Kingdom.

Website: barberboss.co.uk

Email: care@barberboss.co.uk

QR-7116

10

11


12

13

14

Welcome to BarberBoss and congratulations on your new purchase! To ensure your product is fully protected, please activate your 12-month warranty within 90 days by visiting: www.barberboss.co.uk/barberboss-product-warranty. Beyond our warranty, you are covered by your country's legal consumer rights, ensuring complete protection for your purchase.

- 1. Nozzle
- 2. Nozzle Rotation Dial
- 3. Power Button
- 4. Mode Selection Button
- 5. LED Indicator
- 6. USB Charging Port
- 7. Nozzle Release Button
- 8. Water Inlet
- 9. Water Tank
- 10. USB Charging Cable
- 11. Standard Tip X3
- 12. Orthodontic Tip X1
- 13. Periodontal Tip X1
- 14. Tongue Scraper Tip X1

Product Details

Product Name	QR-7116 Travel Water Flosser
Model Number	QR-7116
Charging Input	5V 1.0A
Water Tank Capacity	118 ML
Charging Time	4 Hours
Usage Time	Up to 50 days

Nozzle Tips: Function and Usage

Nozzle Type	Usage	Recommended Replacement Interval
Standard Tip	Ideal for everyday cleaning and removing food debris between teeth and gums.	Every 3 Months
Orthodontic Tip	Designed for cleaning around braces, brackets, and other orthodontic devices.	
Periodontal Tip	Gently cleans deep below the gumline to treat gum disease and improve health.	
Tongue Scraper Tip	Helps remove bacteria and debris from the tongue, promoting fresher breath.	

Charging Instructions:

- Make sure the device is switched OFF before charging.

- Charge the product fully before you use it for the first time and after a long period of disuse. For optimal performance, it is best to use the device when fully charged.
- During Charging, three LED lights flash from Right to Left. The LED Indicator shows charging status:
 - Three Green Lights Solid: Fully Charged
 - Two Green Lights On: Half Power, Device is half charged.
 - Last Light On: Low Battery
 - Last Light Flashes: Needs to Charge
 - Green light blinking: Charging
- Connect provided USB cable to device and plug into a 5V 1A USB adapter or power source (laptop/computer/power bank).

How to Use:

- Insert the tip, then press the power button to begin irrigating; once finished, press the power button to turn off the device.
- The device has a 1-minute auto-stop timer or can be turned off manually by pressing the power button again.
- Pressure Mode Settings: This water flosser offers 3 pressure modes: Gentle (for beginners), Standard (daily cleaning), and Pulse (deep, powerful clean) which can be easily adjusted by pressing the Mode Selection button below LED Indicator.
- To remove the nozzle, press the nozzle tip release button and gently pull the nozzle out.

Precautions:

- Charge only with 5V 1A certified adapter.
- Do not overcharge or charge near water.
- Do not use if damaged or malfunctioning.
- Store between 5°C and 35°C. Avoid damp areas.
- Keep out of reach of unsupervised children under 8 years of age.

Cleaning:

- Power off and disconnect before cleaning.
- Empty the water tank after each use and run the device for 2-3 seconds to drain any remaining water and prevent bacteria buildup.
- Ensure the water tank is completely dry before storing the device for long periods.

Customer Service

- Detailed Instructions, Cleaning, Troubleshooting, Environmental Protection and Recycling: <https://www.barberboss.co.uk/a/FAQ>
- Replacement parts & accessories: www.barberboss.co.uk/shop
- Email: care@barberboss.co.uk
- Phone/WhatsApp: 0800 612 6688 (UK), +44800 612 6688 (International)
- Tutorials: youtube.com/@BarberBossUK

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