*BARBER***BOSS**

QR-7301


Importer: Elite Capital Limited

Address: Studio 4 Clarks Courtyard, 145 Granville Street,
Birmingham, B1 1SB, United Kingdom.


Website: barberboss.co.uk

Email: care@barberboss.co.uk

QR-7301

Brushing Modes

3 Intensity Settings

3 Switch between soft brushing and full power when you need it.

White Mode: Sonic vibration with alternate amplitudes for stain removal and whitening.

Sensitive Mode: Extra gentle mode for sensitive teeth and gums. Recommended mode for beginners.

Gum Care Mode: Long amplitude vibration for gum massage.

Clean Mode: Designed to target and remove plaque for a deep, thorough clean.

Polish Mode: Strong vibration for teeth whitening and polish.

Welcome to BarberBoss and congratulations on your new purchase! To ensure your product is fully protected, please activate your 12-month warranty within 90 days by visiting: www.barberboss.co.uk/barberboss-product-warranty. Beyond our warranty, you are also covered by your country's legal consumer rights, ensuring complete protection for your purchase.

1. Handle
2. Power Button
3. Mode Button
4. Mode Intensity Indicator
5. Mode Light
6. Charging Base
7. Standard Brush Heads x4
8. Mini Brush Heads x2
9. Travel case

Product Details

Product Name	QR-7301 Electric Toothbrush
Model Number	QR-7301
Charging Input	5VDC 0.5A
Charging Time	20 Hours
Usage Time	Up to 180 days

Brush Head: Function and Usage


Head Type	Usage	Recommended Replacement Interval
Standard Head	Medium hardness, DuPont bristles ideal for daily cleaning for most users.	Fading bristles remind you to replace the brush head every 3 months for optimal performance and hygiene.
Mini Head (with silicone protect)	Precise cleaning in hard-to-reach areas, the soft silicone guard offers extra protection for sensitive gums and teeth.	

How to Use:

- Press the Power button to turn on the device.
- Pressure Mode Settings: This toothbrush has 5 pressure modes, easily adjustable by pressing the Power button to change to a new mode. Each mode has 3 intensity levels which can be adjusted by pressing the Mode Selection button below the Power button.
- Turn Off: The toothbrush has a 2-minute smart timer and turns off automatically, or you can also press the power button to stop it anytime.
- Pressure Sensor: Excessive pressure on the bristles triggers rapid vibrations and flashing lights as a warning.

Charging Instructions:

- Make sure the device is switched OFF before charging.
- Fully charge it before first use or after long storage for best performance. To charge, plug the charging cord into a 5V USB adapter and place the handle on the charging base.
- When you turn off the toothbrush, the number of mode lights shows how much battery is left.
- Low Power Indicator: When the “Polish” mode light flashes, it means the battery is low and needs charging soon.

Precautions:

- Charge only with DC 5V certified adapter.
- Do not overcharge or charge near water.
- Do not use if damaged or malfunctioning.
- Store between 5°C and 35°C. Avoid damp areas.
- Keep out of reach of unsupervised children under 8 years of age.

Cleaning:


- Power off and disconnect before cleaning.
- After each use, remove the brush head and rinse it thoroughly. Rinse the handle, especially the brush head connection, at least once a week to maintain cleanliness.

Customer Service


- Detailed Instructions, Cleaning, Troubleshooting, Environmental Protection and Recycling: <https://www.barberboss.co.uk/a/FAQ>
- Replacement parts & accessories: www.barberboss.co.uk/shop
- Email: care@barberboss.co.uk
- Phone/WhatsApp: 0800 612 6688 (UK), +44800 612 6688 (International)
- Tutorials: youtube.com/@BarberBossUK

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